

SPECIFICATIONS & ORDERING INFORMATION

Model No.	Vol. Range	Increment	ISO 8655 Specifications				Microlit Specifications			
			Accuracy		CV		Accuracy		CV	
			±%	± ml	±%	± ml	±%	± ml	±%	± ml
BEAT-2.5	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005	0.6	0.015	0.2	0.005
BEAT-5	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010	0.5	0.025	0.1	0.005
BEAT-10	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020	0.5	0.050	0.1	0.010
BEAT-30	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060	0.5	0.150	0.1	0.030
BEAT-60	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120	0.5	0.300	0.1	0.060
BEAT-100	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200	0.5	0.500	0.1	0.100