

SPECIFICATIONS & ORDERING INFORMATION

Electronic Micropipettes

Model	Volume Range μl	Inc. μl	Acc.		CV	
			$\pm\%$	$\pm\mu\text{l}$	$\pm\%$	$\pm\mu\text{l}$
ME-10	0.2-10	0.05	1.2	0.12	0.8	0.08
ME-120	5-120	0.5	0.8	0.96	0.3	0.36
ME-300	20-300	1.0	0.8	2.4	0.3	0.9
ME-1000	50-1000	5.0	0.8	8	0.3	3
ME-5000	100-5000	10.0	0.8	40	0.3	15

Error limits according to the nominal capacity (=maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20° C, and with smooth, steady operation. The Error limits are well within the limits of DIN EN ISO 8655.2.