

SPECIFICATIONS & ORDERING INFORMATION

Miniature Fixed Volume

Model	Vol. μl	Acc. \pm%	CV \bullet%
MM-5	5	1.5	1.0
MM-10	10	1.0	1.0
MM-20	20	0.5	0.5
MM-25	25	0.5	0.5
MM-50	50	0.4	0.3
MM-100	100	0.3	0.3
MM-200	200	0.5	0.5
MM-250	250	0.5	0.5
MM-500	500	0.5	0.5
MM-1000	1000	0.5	0.5