

Fully Autoclavable Fixed Volume Micropipettes

Model	Vol. μl	A		CV		Model	Vol. μl	A		CV	
		$\pm\%$	$\pm\mu\text{l}$	$\pm\%$	$\pm\mu\text{l}$			$\pm\%$	$\pm\mu\text{l}$	$\pm\%$	$\pm\mu\text{l}$
RBO-F-1	1.0	5	0.05	5	0.05	RBO-F-200	200.0	0.6	1.2	0.2	0.4
RBO-F-2	2.0	4	0.08	2	0.04	RBO-F-250	250.0	0.6	1.5	0.2	0.5
RBO-F-5	5.0	2	0.1	1	0.05	RBO-F-500	500.0	0.6	3	0.2	1
RBO-F-10	10.0	1	0.1	0.5	0.05	RBO-F-1000	1000.0	0.6	6	0.2	2
RBO-F-20	20.0	0.8	0.16	0.4	0.08	RBO-F-2000	2000.0	0.6	12	0.2	4
RBO-F-25	25.0	0.8	0.2	0.4	0.1	RBO-F-5000	5000.0	0.6	30	0.2	10
RBO-F-50	50.0	0.8	0.4	0.4	0.2	RBO-F-10000	10000.0	0.4	40	0.2	20
RBO-F-100	100.0	0.6	0.6	0.2	0.2						