

Fully Autoclavable Variable Volume Micropipettes

Model	Volume Range	Inc. μl	A		CV		Model	Volume Range	Inc. μl	A		CV	
			$\pm\%$	$\pm\mu\text{l}$	$\pm\%$	$\pm\mu\text{l}$				$\pm\%$	$\pm\mu\text{l}$	$\pm\%$	$\pm\mu\text{l}$
RBO-2	0.2 - 2.0 μl	0.002	2	0.04	1.2	0.024	RBO-200	20 - 200 μl	0.2	0.6	1.2	0.2	0.4
RBO-10	0.5 - 10 μl	0.02	1	0.1	0.5	0.05	RBO-1000	100 - 1000 μl	1.0	0.6	6	0.2	2
RBO-20	2 - 20 μl	0.02	0.8	0.16	0.4	0.08	RBO-5000	0.5 - 5 ml	10.0	0.6	30	0.2	10
RBO-50	5 - 50 μl	0.1	0.8	0.4	0.4	0.2	RBO-10000	1 - 10 ml	20.0	0.6	60	0.2	20
RBO-100	10 - 100 μl	0.2	0.6	0.6	0.2	0.2							