

**SPECIFICATIONS &
ORDERING INFORMATION****ISO 8655 Specifications****Model
No.****Vol.
Range****Increment****Accuracy****±% ± ml****CV****±% ± ml**

BTDR-1

0.25-2.5 ml

0.05 ml

0.6

0.015

0.2

0.005

BTDR-2

0.5-5 ml

0.1 ml

0.6

0.030

0.2

0.010

BTDR-3

1-10 ml

0.2 ml

0.6

0.060

0.2

0.020

BTDR-4

2.5-30 ml

0.5 ml

0.6

0.180

0.2

0.060

BTDR-5

5-60 ml

1.0 ml

0.6

0.360

0.2

0.120

BTDR-6

10-100 ml

2.0 ml

0.6

0.600

0.2

0.200

BTDR-7

50-400 ml

5.0 ml

0.6

2.400

0.2

0.800