

# SPECIFICATIONS & ORDERING INFORMATION

Model No.	Vol. Range	Increment	Accuracy		CV	
			±%	± ml	±%	± ml
ULT-2.5	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005
ULT-5	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010
ULT-10	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020
ULT-30	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060
ULT-60	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120
ULT-100	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200

Error limits (Accuracy & Coefficient of variation) according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth, steady operation. The error limits are well within the limits of DIN EN ISO 8655-2.